

**Technical datas**

- Weight** 3,1 kg with battery.
- Size** Height 27 cm;  
Width 10 cm;  
Length 24 cm.
- Oxygen Flow** Pulse dose delivery system. Four flow settings from 1-6.
- Power Sensing** Universal Power Supply with autosensing for worldwide use.  
In AC mode for 100-240V; in DC mode for 12V use in Car/airplane/boat/camper.
- Battery** Duration 3.5 - 4 hours\*. 7 - 8 hour battery available as an option.
- Noise Level** Very low at 38 Decibels.\*
- Warranty** 2 year limited warranty.\*\*
- Operation** Simple control functions; clear LCD display.
- Use** Stationary, portable and travel. Clinically effective for day and night use.

**Already included**



**Lithium Ion Battery**  
The battery lasts for up to four hours. When it's plugged in (AC or DC), the battery automatically recharges. Eight hour batteries are also available.



**Universal Power Supply**  
When plugged into an AC or DC power outlet, the INOGEN ONE G2 delivers an unlimited supply of oxygen.




**Carry bag**  
Designed specifically for the INOGEN ONE G2. Take advantage of your new independence and mobility.



**Nasal Cannula**  
A device for dispensing oxygen (length 2.1 m).

\* At setting 2, in typical home use.  
\*\* See limited warranty statement for complete details.

# INOGEN ONE G2



**Helgas is the official INOGEN distributor**  
**FOR MORE INFORMATIONS: +39 0161.1896507**

HELGAS di Grazia Presti - via Facchinetti 3/3A, 13100, Vercelli (IT)  
Tel +39 0161.1896507 - Fax +39 011.19837151 - E-mail info@oxypoint.com - Web www.oxypoint.com

**Oxygen: wherever, whenever.**

# Restore your freedom with INOGEN ONE G2

Oxygen therapy should not have to limit you in any way. INOGEN ONE G2 offers just that, no worries about running out of oxygen, no need to keep an eye on the clock, no need to carry heavy equipment around. Simply use the INOGEN ONE G2 wherever you are, wherever you go, whenever you want.

## Wherever you are

It sits discreetly at your feet, behind a chair, or next to the bed and it's light enough to carry over your shoulder. All of your oxygen needs can be met with this one small device, day or night, home or away. When using the INOGEN ONE G2 you will enjoy all the benefits of this compact and lightweight device, producing just a very light sound, hardly noticeable to you or those around you. And, it's easy to operate too, just a couple of buttons. Embrace the simplicity of the INOGEN ONE G2 and see how it can make your life easier.

## Wherever you go

Restoring your freedom is definitely a large part that makes the INOGEN ONE G2 so very special. So go ahead, jump in the car to run errands or take a weekend trip to see your family, or make that long journey you'd never thought you'd be able to make. Head home when you want to, not when a limited oxygen supply forces you to.

Most airlines allow you to use the INOGEN ONE G2 when you fly, making it easier than ever to reclaim your independence.



## Whenever you want

INOGEN ONE G2 is built to withstand long term and heavy duty use.

Expect it to run for years on end and it will need full maintenance only after it has been running for 30,000 hrs.

Your INOGEN ONE G2 oxygen system can be delivered to you quick and easy, anywhere in Europe, through [www.oxypoint.com](http://www.oxypoint.com). Please contact us!



## Extra accessories



Back pack



Stand alone charger



Cable for airplane seat



8/16 cells additional battery



Particulate filter



Antibacterial filter

INOGEN ONE G2 is approved by every airline company.

