

INOGEN ONE G2



**12 CELLS
BATTERY**



**24 CELLS
BATTERY**

WEIGHT	3,1 KG
DIMENSIONS	27,4 × 10,1 × 24,1 cm
NOISE	Less than 38 dBA (as packaged) on setting 2
12 CELLS BATTERY	Up to 4 hours on setting 2
24 CELLS BATTERY	Up to 8 hours on setting 2
OXYGEN CONCENTRATION	90% - 3% /+ 6% at all settings
NIGHT FUNCTION	AUTOMATIC
ACTIVATION PRESSURE	12 cm H ² O
FLOW CONTROL SETTINGS	6 settings: 1 to 6
ENVIRONMENTAL RANGES INTENDED FOR USE	Temperature: 41 to 104°F (5 to 40°C) Humidity: 0% to 95%, non-condensing Altitude: 0 to 10,000 ft (0 to 3048 meters)
ENVIRONMENTAL RANGES INTENDED FOR SHIPPING AND STORAGE	Temperature: -13 to 158°F (-25 to 70°C) Humidity: 0% to 95%, non-condensing Store in a dry environment Altitude: 0 to 10,000 ft (0 to 3048 meters)
DATA LOG	1 year / USB o RF
EC CERTIFICATE	CE0123
MEDICAL DEVICE NUMBER	473530

INOGEN G2		BOLO, ML			
SETTING	FLOW ml/min	10 BPM	15 BPM	20 BPM	25 BPM
1	180	18	12	9	7
2	360	36	24	18	14
3	540	54	36	27	22
4	720	72	48	36	29
5	900	90	60	45	36
6	1080	108	72	54	43

BPM: Breaths Per Minute